Dr. Fog

Aerosol **digitale**

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PARTS

- | | |
|---------------------------------|-------------------------------|
| 1. LED screen display | 10. Mouthpiece |
| 2. Transformer plug | 11. Nebulizer |
| 3. On/off button | 11a. Nebulizer cover |
| 4. Compressor (main unit): | 11b. Nebulizer top section |
| 5. Nebulizer support | 11c. Nebulizer bottom section |
| 6. Air filter and cover housing | 12. Air tube |
| 7. Air tube fitting | 13. Air filters |
| 8. Adult mask | 14. Transformer |
| 9. Child mask | 15. Carry bag |

SUPPLIED WITH

Compressor (main unit); nebulizer; adult mask; child mask; mouthpiece; air tube; five air filters; adaptor; instruction manual; carry bag.

If a part should be missing, please contact the retailer.

NOTE: Parts in contact with the skin are tested according to ISO 10993 standards and do not have present any allergy risk.

NOTE: Accessories (nebulizer, adult mask, child mask, mouthpiece, air tube) are for re-use by a single patient; they are not to be shared with other users.

Thank you for choosing the Dr. Fog digital aerosol from J BIMBIR, designed for use with children.

The device is an accessory designed to be connected with a nebulizer so as to constitute a nebulization system that delivers prescribed medicines to patients affected by cystic fibrosis, asthma, BPCO, and other illness or problems linked to the respiratory system. This device must be used under the supervision of a doctor. This manual should be read with care to familiarise with the characteristics of the product. It is necessary always to avoid any use of the product other than for its intended purpose.

OPERATING PRINCIPLE

The compressed air aerosol is a treatment device that uses the Venturi principle to nebulize medicines using compressed air. The compressed air is conveyed through a converging nozzle, where it accelerates and comes out at high speed, creating a vacuum (Venturi effect). The vacuum extracts the liquid from a tank through a cylindrical channel and it conveys it in the outgoing air flow, formed by a nozzle, where it is mixed with air and impacts against a rigid surface. This process uses the energy of the air flow to convert the liquid into small droplets, known as aerosol. Once they reach the user, these nebulized liquid particles are able to penetrate the lungs effectively (Fig. E).

According to medical research, the size of the particles plays an important role in the way they deposit in the lungs, together with particle speed and sedimentation time.

As can be seen in the figure, larger particles, from 10 to 15 μm , mainly deposit in the upper airways, while particles between 5 and 10 μm in size reach the primary bronchi, while those that are 1–5 μm penetrate the lower airways and the peripheral lungs.

The advantage of inhaled aerosol treatment is that compared to administration by mouth, inhaled medicines their action is faster (Fig. F). (*Taken from "A GUIDE TO AEROSOL DELIVERY DEVICES FOR RESPIRATORY THERAPISTS" - 5th edition*).

HOW TO USE

This device has been designed to produce compressed air needed for the operation of a nebulizer to produce therapeutic aerosol particles for the treatment of breathing problems.

RECIPIENTS

The device is intended for use by medical personnel (such as doctors, nurses, therapists) and adult patients.

The user must be able to understand the general operation of the device and the contents of the instruction manual.

PATIENTS


The device is intended for adult and elderly patients; newborns and children can only use the device under adult supervision.

PERFORMANCE INDICATOR

Table 1: the test has been performed with a 0.9% saline solution at a temperature of 25 ± 3 °C and relative humidity of $60 \% \pm 5 \%$. Equivalent volume of nebulized particles.

Performance indicator table

Nebulization capacity of the nebulizer	$\geq 0.3\text{ml/min}$
Medicinal fluid temperature	$\leq 40\text{ }^{\circ}\text{C}$
Noise emissions during operation	$\leq 65\text{dB}$
Diameter of nebulized particles	MMAD: $3.8\mu\text{m} \pm 25\%$
Residual quantity	$< 1\text{mL}$

NOTE: the horizontal line indicates the diameter of the particles; the data represent the distribution logarithm. The vertical line shows the cumulative percentage in volume, corresponding to the upward trend of the curve; the vertical line on the right indicates a percentage interval in volume terms, corresponding to the histogram or mobile curve.

SAFETY INFORMATION

NOTE: READ ALL INSTRUCTIONS CAREFULLY BEFORE USE.

The basic precautions listed below are to be followed when using an electrical product.

Important: failure to read and abide by all warnings and precautions may lead to personal injury or damage to the device.

CONTRAINDICATIONS

- To be used exclusively on conscious patients.
- Use is prohibited on patients with seriously low levels of oxygen or severe respiratory distress.
- In case of chronic illness, hospitalisation or intensive care stays in the previous six months, the device may only be used under medical supervision.
- This device is not suitable for use in anaesthesia breathing systems or in pulmonary ventilation systems.
- This device cannot be used with pentamidine.

SIDE EFFECTS

- Systemic administration of aminoglycosides, including AMK, is at times associated with serious side effects, including hearing toxicity, i.e., ototoxicity, vestibular toxicity and renal toxicity.
- Oral discomfort.
- Dysphonia.

WARNINGS AND PRECAUTIONS

Warnings

- Pentamidine is not an approved medicine for use with this device! This device cannot be used in anaesthesia breathing systems or in pulmonary ventilation systems.

- Not suitable for use in the presence of flammable anaesthetic mixtures with air, oxygen or nitrous oxide.
- Only use medicines prescribed by your doctor and follow the instructions concerning dose, duration and frequency of treatment.
- Children must always use this device under adult supervision.
- The device contains small parts that are easily swallowed.
- Be aware of the risk of strangulation from cables and flexible tubes.
- Never place the device in such a way that it is difficult to unplug it from the power supply.
- The main compressor unit is re-usable and can be used by more than one person after suitable cleaning.
- The nebulizer accessories are for use on a single patient only.
- This product is designed for use with medicines in solution form.
- Never use with medicines containing fats or suspended particles, including herb-based medicines. Common medicines for delivery by nebulizer are recommended, as instructed by your doctor.
- The device nebulizes for 10-20-30 minutes per session, using a quantity of medicine from 2 to 8 ml. However it can be used as instructed by your doctor.
- As far as concerns the type, dose and treatment regime of the medicines concerned, always follow the doctor's advice.

Important: unauthorised use of medicines via aerosol inhalation may put your life at risk!

- This device is not intended as life support or to maintain vital functions.
- Never share the accessories with other people. Allowing more than one person to use the device risks spreading infectious diseases.
- This nebulizer is exclusively for use in respiratory therapy; any other use of the device is improper and therefore, dangerous.
- Please use the accessories (nebulizer, adult mask, child mask, mouthpiece, etc.) supplied by the manufacturer. The accessories have been tested compliant with the ISO 10993-1 standard.

Product Warnings

- To avoid electric shocks:
 - Keep the device away from water.
 - Never place the power cord or device in liquid.
 - Never use the device while bathing.
 - Never try to recover a device that has fallen into water: disconnect it from the power socket immediately.
- Never use the device if any parts are damaged (including power cord), if it has been dropped or if it has been immersed in water. Contact your assistance service immediately for an inspection and if necessary, repairs.
- The device must not be used in the presence of fumes or flammable gases.
- Keep the air intakes free from obstruction. Never place the device on a soft surface that could block the openings.
- In case of anomalies, cease all use of the device until it has been checked and repaired.
- The device must never be left unattended when connected to the power socket.
- Never angle or shake the appliance while it is operating.
- Disconnect the device from the power socket before cleaning and after every use.
- The use of solutions, suspensions or emulsions other than those recommended by the manufacturer, and in particular a suspensions and/or solutions with high viscosity, may alter the granulometric delivery curve, the mass median aerodynamic diameter (MMAD), and the aerosol production and/or capacity, which could all then differ from the values indicated by the manufacturer.

Precautions for use

- Connect the device to a power socket with the correct voltage for the model being used.
- Never leave the device to operate unsupervised.
- Never use the device if the power cord or plug are damaged, if the device has fallen in water, or if it is not operating correctly. Contact an assistance service to carry out the repairs.
- In case of anomalies, cease all use immediately and until the device has been checked

and repaired.

- Always disconnect the device from the power socket immediately after use.
- Never obstruct the air intakes of the compressor or use it in places where the air intakes could be obstructed.

Storage Precautions

- Never store the device where it is exposed to direct sunlight, high temperatures or elevated moisture levels.
- Keep the product out of the reach of small children.
- Always keep the device unplugged from the power socket when not in use.
- The time needed for the compressor to cool down from maximum storage temperature between uses until it is ready for its intended use is 20 minutes.

Cleaning Precautions

- Never immerse the device in water. This could damage it.
- Unplug the device from the power socket before cleaning.
- Clean all parts that require it after every use, following the instructions in this instruction manual.

USING THE DEVICE

Before use, make sure that the mask, mouthpiece and air filter are all clean.

Power Supply:

Connect the adaptor to the connector on the compressor and to the power socket. Never perform this operation with wet hands.

Insert the medicine (Fig. A)

Unscrew the top part of the nebulizer and lift it off.

Important: take care that the nebulizer cone does not fall onto the floor. The nebulizer

cone is secured to the top part of the nebulizer by two hooks and can be removed using a certain amount of force.

Pour the correct amount of medicine, as prescribed, into the bottom section of the nebulizer.

Nebulizer Assembly

Turn the top section of the nebulizer clockwise to screw it onto the bottom section until the nebulizer is tightly closed (Fig. A).

Fit the required breathing accessory: mouthpiece (Fig. B) or mask (Fig. C).

Fitting the Air Tube

Fit one end of the air tube to the air tube connector on the compressor, pushing and firmly turning it into place (Fig. D).

Connect the other end of the air tube to the nebulizer, firmly pushing it into the connector on the bottom of the nebulizer (Fig. D).

Temporarily use the nebulizer support (5) if fitted to the compressor to perform this operation.

Important: keep the nebulizer upright while fitting the air tube. Failure to do so could cause the medicine to leak from the nebulizer.

Important: never use the device if the air tube is bent.

Using the device

Hold the nebulizer in one hand, as shown in Figure B.

Important: never hold the nebulizer at an angle above 45°. If this is not the case, the medicine could pour into the mouth.

Press the on/off button (3) to start the device.

The LED screen (1) will now read “30 minutes” and nebulizing will start.

Press the on/off button again to select 20 minutes, or press once more to select 10 minutes.

Press the button again to switch off the device.

The device will operate for the selected time or for as long as there is medicine to dispense,

and the remaining minutes will be displayed as a countdown.

10/20/30 minute timer function

Thanks to the timer function, it is possible to set the treatment time.

1. 8 ml of medicine is delivered in approx. 30 ± 5 minutes.
2. 6 ml of medicine is delivered in approx. 20 ± 3 minutes.
3. 3 ml of medicine is delivered in approx. 10 ± 2 minutes.

DELIVERING TREATMENT

The medicine is to be inhaled as per doctor's instructions.

The difference in the inhalation accessories are as follows.

- **Using the mask:** bring the mask up to the face, covering nose and mouth (Fig. C). Adjust the elastics to secure it comfortably. Breathe in and out through the mask.
- **Using the mouthpiece:** bring the mouthpiece up to the mouth and gently close the lips around it to allow the regular flow of the nebulized medicine. Breathe in and out through the mouthpiece (Fig. B).

End of inhalation:

- Disconnect the air tube from the nebulizer: hold the air tube connector and gently pull it downwards. Check the air tube and make sure it contains no traces of condensation or moisture inside the tube.
- Follow the instructions below:
 - Make sure that the other end of the air tube is still connected to the compressor.
 - Switch on the device. The compressed air will pass through the tube and expel all traces of moisture or condensation.
 - Switch off the device.
- Disconnect the air tube from the compressor: take the air tube connector and take it out of the air tube connector on the compressor.
- Disconnect the adaptor from the power socket and compressor.

CLEANING AND DISINFECTION

CLEANING

Be sure to clean and store the device after every use.

- Remove the accessories from the nebulizer.
- Disconnect the air tube from the nebulizer.
- Gently unscrew the top part of the nebulizer.
- Empty out any remaining medicine.
- Carefully wash the parts in water.
- Dry by hand or air dry in a clean environment, using a soft, clean cloth.
- Refit the nebulizer and place it in the carry bag.
- Clean the main unit with a soft cloth and a little water or mild detergent. Never use abrasive detergents.

Make sure that the internal parts of the device do not come into contact liquids and that the power plug is disconnected.

DISINFECTION

Disinfect the mask or mouthpiece after every treatment during the day.

Use a normal disinfectant, commonly found on sale, following the manufacturer's instructions; this procedure is also recommended to disinfect masks and mouthpiece.

Immerse the parts in the disinfectant solution for the instructed time.

Remove the parts and dispose of the solution.

Rinse the parts in hot tap water; shake off any excess water and dry in a clean environment.

Never use a microwave oven to dry the parts.

Never use autoclaves, EOG gas disinfection systems, or low-temperature plasma sterilizers to disinfect the device.

MAINTENANCE

No institute or individual is authorised to perform maintenance or repairs on the product. In case of doubt or questions about the device, please contact the manufacturer or retailer for more information.

Users must never try to repair the device or its accessories. Always contact the retailer for any repair requirements.

Unauthorised persons are not permitted to open the device and any such opening will void any rights to the warranty.

Clean/disinfect the device after every use: for more details about the method, see the 'Cleaning and Disinfection' section.

The compressor and accessories must not undergo maintenance or repairs while being used by a patient.

WARNING: Changes to this device are not permitted!

AIR FILTER REPLACEMENT

- Never use the device without an air filter installed.
- Never use cotton or other materials.
- Never wash or clean the air filter.
- Replace the filter every 30 days or when it begins to turn grey. Only use filters provided by your retailer.
- How to replace the filter:
 - Remove the filter housing cover.
 - Replace the used filter with a new one.
 - Refit the filter housing cover.

Important:

- Do not attempt to wash the filter for further use.
- The air filter must not undergo maintenance or repair while being used by a patient. Only use original filters! Do not use the device without a filter installed.

TROUBLESHOOTING

Problem	Possible cause	Solution
No nebulization or low nebulization speed when the device is on.	No medicine in the nebulizer. Too much or too little medicine in the nebulizer.	Insert the correct amount of the prescribed medicine into the nebulizer.
	The nebulizer is incorrectly assembled.	Make sure the nebulizer is correctly assembled and ensure the accessories are correctly connected.
	The nebulizer cone is not secured to the top part of the nebulizer or it is incorrectly positioned.	Make sure that the nebulizer cone is correctly secured to the top part of the nebulizer.
	The nebulizer cone nozzle is blocked.	Clean and disinfect the nebulizer to remove the obstruction.
	The nebulizer is being held at an incorrect angle.	Hold the nebulizer correctly. Never hold the nebulizer at an angle above 45°.
	The air tube is incorrectly connected.	Make sure that the air tube is correctly connected to both the compressor and the nebulizer.
	The air tube is bent, damaged or blocked.	Make sure the air tube is not bent, twisted or squashed. Make sure that the air tube is undamaged. If damaged, the tube must be replaced.
	The air filter is dirty.	Replace the air filter with a new, clean one.

Problem	Possible cause	Solution
The device does not switch on when the on/off button is pressed.	The transformer is not connected to the compressor.	Connect the transformer to the compressor and press the on/off button.
The device overheats.	The compressor is covered.	Never cover the compressor with anything.
	The compressor has been working continuously for more than 30 minutes.	Limit use to 30 minutes at a time and leave the device for 20 minutes before using it again.
The device is unusually noisy.	The air filter is not correctly inserted.	Make sure that the filter cover has been correctly fitted and that it is not blocked.

Note: if the recommended solution does not solve the problem, do not try to repair the device: no part of the device can be repaired by the user. Take the device to an authorised retailer or distributor.

TECHNICAL SPECIFICATIONS

Product name	Digital Aerosol
Model	Code: 12634 (CDC300S)
Accessories	Nebulizer, Air tube, Adult mask, Child mask, Mouthpiece, 5 air filters, transformer, carry bag, instruction manual
Power Supply	AC 100v – 240V, 50Hz/60Hz
Consumption	Power supply: Input: 100-240 V, 50/60 Hz, 0.35 A max. Output: 12 V DC, 1°
Electric shock protection	Class II device
Electric shock protection grade	BF type part applied
Operating method	Continuous operation
Nebulization speed	>0,3 ml/min
Maximum nebulizer capacity	8 ml
Sound level	< 65dB(a)
Particle size	MMAD<3.8um um±25%
Compressor pressure interval	8.0-19.0 psi
Working pressure interval	5.8-10 psi
Capacity	≥3.5 l/min
Operating temperature	5-40°C
Operating humidity	15-80% HR
Operating air pressure	70-106 KPa
Storage temperature	-20-60°C
Storage humidity	10-90% HR
Storage air pressure	70-106 KPa
Dimensions	138x103x64mm
Weight	Approx. 300 g (without nebulizer or air tube)
Protection against harmful water or particulate damage	IP21
Useful life	5 years

WARRANTY

This device has a 2-year warranty from the date of purchase. The warranty is valid only with the original purchase receipt.

- The 3-year warranty is valid only for the compressor.
- The accessories such as nebulizer, masks, mouthpiece, air tube and filters have a warranty of one year against defects in materials and workmanship, but not against deterioration caused by normal wear and tear. In the event that replacements are needed, we recommend using only original accessories.
- Opening or tampering with the device voids the warranty.
- The warranty does not cover damage caused by improper use, accidents, falls, immersion in liquids or failure to abide by the instructions for use.

NOTE: If the device needs any repairs, please contact your local distributor. In the meantime, please keep the packaging (box) in a safe place so that it can be used for shipping.

DISPOSAL

Dispose of the device in compliance with current regulations in the place of use. In EU countries, dispose of the device at a public waste disposal point – WEEE Directive 2002/96/EC. In case of any doubts, please contact your local authority for information on waste disposal.

ELECTROMAGNETIC COMPATIBILITY

1st WARNING: use of this equipment is to be avoided close to or stacked on other equipment, since this could lead to malfunctions. If such use is required, keep this and the other equipment under supervision to ensure they are operating correctly.

2nd WARNING: the use of accessories, transformers or cables other than those specified or supplied by the manufacture of this appliance may cause an increase in electromagnetic emissions or a reduction in electromagnetic immunity of the appliance, causing incorrect operation.

3rd WARNING: portable and mobile radio-frequency communication equipment (including peripherals such as aerial cables and external aerials) cannot be used at a distance of less than 30 cm from any part of the device, including the cables specified by the manufacturer. If this is not the case, it may lead to deterioration in equipment performance.

Table 1 • Declaration – Electromagnetic emissions

Emissions test	Conformity
RF emissions (CISPR 11-3-2)	Group 1
RF emissions (CISPR 11-3-2)	Class B
Harmonic emissions (IEC 61000-3-2)	Class A
Voltage fluctuations/flicker (IEC 61000-3-3)	Compliant

Table 2 • Declaration – Electromagnetic immunity

Immunity test	Test level IEC 60601	Conformity level
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air
Electrical fast transients/burst IEC 61000-4-4	±2kV for power supply line	±2kV for power supply line
Peak (IEC 61000-4-5)	±0.5kV, ±1 kV line to lines	±0.5kV, ±1 kV line to lines
Voltage dips, short interruptions and voltage variations in the power network IEC 61000-4-11	0% UT, 0.5 cycles, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT, 1 cycle and 70% UT, 25/30 mono-phase cycles at 0° 0% UT, 250/300 cycles	0% UT, 0.5 cycles, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT, 1 cycle and 70% UT, 25/30 mono-phase cycles at 0° 0% UT, 250/300 cycles
Industrial frequency magnetic field (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m

NOTE: UT is the mains voltage in AC before application of the test level.

Table 3 • Declaration - Electromagnetic immunity

Immunity test	Test level IEC 60601	Conformity level
Conducted RF IEC 61000-4-6	3V 0.15 MHz at 80 MHz 6 V ISM and amateur radio bands between 0.15 MHz and 80 MHz	3V 0.15 MHz at 80 MHz 6 V ISM and amateur radio bands between 0.15 MHz and 80 MHz
Radiated RF IEC 61000-4-3	10V/m 80 MHz at 2.7 GHz	10V/m

**Table 4 • Declaration – IMMUNITY to proximity effects
from RF wireless communication equipment**

Immunity test	Test level IEC 60601				Conformity level
	Test frequency	Modulation	Maximum power	Immunity level	
Radiated RFIEC 61000-43	385 MHz	**Pulse modulation: 18Hz	1.8W	27 V/m	27 V/m
	450 MHz	*FM + 5Hz deviation: sinusoid at 1kHz	2 W	28 V/m	28 V/m
	710 MHz 745 MHz 780 MHz	** Pulse modulation: 217Hz	0.2 W	9 V/m	9 V/m
	810 MHz 870 MHz 930 MHz	** Pulse modulation: 18Hz	2 W	28 V/m	28 V/m
	1720 MHz 1845 MHz 1970 MHz	** Pulse modulation: 217Hz	2 W	28 V/m	28 V/m
	2450 MHz	** Pulse modulation: 217Hz	2 W	28 V/m	28 V/m
	5240 MHz 5500 MHz 5785 MHz	** Pulse modulation: 217Hz	0.2 W	9 V/m	9 V/m

NOTE * – As an alternative to FM modulation, it is possible to use a pulse modulation of 50% at 18 Hz, since although it does not represent real modulation, it would constitute a worst-case scenario.

NOTE ** – The carrier must be modulated using a square wave signal with 50% work cycle.

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Shenzhen Fitconn Technology Co., Ltd
Wellkang Ltd, (www.CE-marking.eu)